Neural implant insertion system using ultrasonic micro-vibration to improve #700.26 penetrating microelectrode array insertion mechanics: In vitro and in vivo evaluations



Innovative motion + Positive outcomes

Introduction

Penetrating electrode arrays provide localized access to neural stimulation and recording sources, and they are a valuable neuroscience research tool. Neuroscientists and researchers developing future clinical applications continually seek to increase information transfer potential by implanting electrode arrays with higher numbers of shanks at greater density. Successful insertion of penetrating electrode arrays into brain tissue remains a significant challenge, and outcomes are often heavily reliant on surgical skill and technique. During insertion, electrode arrays apply forces to the neural tissue resulting in significant compression (dimpling) at the implant site and strain on local tissues. Shallow insertions (< 1 mm) targeting the upper cortical layers are especially difficult. Electrode arrays may need to be over-inserted past the target depth to ensure puncture of meningeal membranes, or be fired in at very high velocity, limiting accuracy and possibly increasing implantation trauma.

NeuralGlider-Cortical[™] uses ultrasonic micro-vibration to reduce the insertion force and dimpling during microelectrode array implantation.

NeuralGlider-CorticalTM System

- Ultrasonic actuator produces axially-directed micro-vibrations in shanks as array is inserted
- Stereotaxis mountable closed loop micro-positioner stage, 0.5 µm resolution
- Omnetics-compatible, detachable coupler for reliable transmission of ultrasonic vibration to the neural implant
- RoHS compliant control box
- Low profile USB microscope camera, with user-positioned mount to the stereotaxic instrument
- LabVIEW-based GUI for control of the micropositioner and actuator while recording position data

Prior Research

- Initial *in vivo* study:
- Utilized custom-built 8-channel tungsten microwire arrays (50 µm diameter with 250 µm spacing)
- Bi-lateral barrel cortex implants
- All insertions through pia (dura removed) at 50 µm/s
- 2-week survival (N=10); recordings and histology
- Ultrasonic micro-vibration of penetrating electrode arrays during insertion substantially reduces insertion forces and dimpling of brain tissue
- Initial short-term studies supported safety of approach; no significant difference in recording quality or histology



nsertion Depth: 850 um 8-ch arrays inserted with (left) and without (right) ultrasonic vibration into rat barrel cortex. Dashed lines: brain surface at initial contact (Blue) and at target insertion depth (Green).

Ryan S. Clement, Ph.D*., Natasha N. Tirko, Ph.D., Jenna K. Greaser, Casey A. Scruggs, Maureen L. Mulvihill, Ph.D., and Roger B. Bagwell, Ph.D.

Actuated Medical, Inc., 310 Rolling Ridge Dr., Bellefonte, PA 16823 *ryan.clement@actuatedmedical.com







Maximum dimple, estimated from analysis of insertion videos. (* p<0.001; n=10)

Improving Insertion Mechanics of Commercial Microwire Arrays



Above: Examples of agar insertion force vs. insertion distance profiles for both Tucker-Davis and MicroProbes arrays

Right: Summary analysis of initial puncture force (1st peak in force vs. insertion distance plot) vs. actuator power.

Insertion Force and Dimpling Reduction: Ex Vivo Evaluations





Distance (mm)

Examples of *ex vivo* insertion force testing sample preparations. Left: sample with large craniotomy on bed of agar for pia-only insertions. Right: Example of preparation with small craniotomies for testing dura insertion approach.

Demonstrated Safety and Efficacy of Ultrasonic Vibration-aided Insertion, Including Through Dura



Video frames corresponding to 250 µm intervals during insertion process were extracted and analyzed. Average brain surface dimple determined from averaging exposed lengths of visible wires relative to tip positions when brain first contacted (Image-J).

* Max dimple estimated by total distance array was advanced as blood obscured video and/or rupture (and subsequent brain relaxation) was not clear.



Distance at Vibration Onset re: Surface (µm)

Tucker-Davis

Substantially more "off-axis" vibration was observed with MicroProbes (top) arrays when actuated in air as compared to Tucker-Davis. Delaying vibration until tips make initial contact avoids off-axis vibration (middle).

Delaying vibration until the tips were partially advanced into the model did not greatly mpact performance.



Ex vivo insertion force summary for insertion of Tucker-Davis arrays through pia only vs. through intact dura (control = non-vibrated). With ultrasonic vibration, insertion through dura generated comparable forces as insertion through pia only. Insertion through small craniotomy with intact dura yielded less dimpling.



Compatibility with Many Array Types

Blackrock

Poly-ethylene glycol (PEG) was used to attach unwired 4x4 array to actuator coupler. Below: Plots of insertion force and dimpling vs. distance for representative vibrated and non-vibrated insertions. *Right:* Video frames from same insertions extracted after advancing 1000 µm. Dashed line marks surface of agar model at initial contact (blue), and at maximum dimple (green).



NeuroNexus

Several NeuroNexus probe configurations have been evaluated with the NeuralGlider system, including both acute and chronic varieties. Reductions in force and dimpling (evaluated in agar model) are 60-80%. Chronic (CM) Style Probe Acute Style Probe





Additional System Capability

Neural Recording During Insertion



With custom angle adapter to couple headstage amp to array, it is possible to obtain neural recordings at points during the insertion process

Conclusions and Future Directions

Acknowledgments

This work was sponsored by the Defense Advanced Research Projects Agency (DARPA) Biological Technologies Office (BTO) Electrical Prescriptions (ElectRx) program under the auspices of Dr. Doug Weber and Dr. Eric Van Gieson (DARPA Contracts Management Office Contract No. HR0011-16-C-0094), as well as the NIH BRAIN Initiative (1R43NS105500-01A1). The content is solely the responsibility of the authors and does not necessarily represent the official views of the Defense Advanced Research Projects Agency or the NIH. Pat. actuatedmedical.com/ip

Neural Implant Insertion System

CORTICAL

	350	
e Ne	300	$\widehat{}$
1	250	un)
	200	ole
	150	jm
		<u>e</u>
、	50	erag
1	0 1.25	Ave



► N-Form

PEG was also used to attach a 2x2 array with portion of ribbon cable to actuator coupler. Insertion force measurement setup with agar test sample resting on force sensor is shown in lower left photo. Insertion force and dimpling could be reduced by 75% and 60%, respectively. PEG can be dissolved with water to release array.







Distance (mm)

Compatibility with TDT's ZIF Connector



Demonstrated compatibility of NeuralGlider with **Tucker-Davis** Technology's (TD1 patented ZIF connector arrays popular among researchers





Ex vivo pilot study demonstrated that NeuralGlider may improve Inscopix GRIN lens insertions (deeper penetration with less brain compression)

+ Ultrasonic micro-vibration by the NeuralGlider neural implant system safely improves insertion of microwire arrays and many other array types on the market 60-80 % reductions in tissue dimpling and insertion force are typical The NeuralGlider system enabled the insertion of microwire arrays through dura, which may simplify future neural implant procedures for some devices - The system is shown to be compatible with many electrode array types; also could have application in the placement of other non-electrode devices + Future *in vivo* studies required for evaluating potential for improving chronic results and for evaluating performance with non-microwire devices

